


ba12217e

# inoLab pH Level 1



Laboratory pH meter

---

**Accuracy when  
going to press**

The use of advanced technology and the high quality standard of our instruments are the result of continuous development. This may result in differences between this operating manual and your instrument.

We cannot guarantee that there are absolutely no errors in this manual. We are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions.

**Warranty  
declaration**

The designated instrument is covered by a warranty of three years from the date of purchase.

The instrument warranty extends to manufacturing faults that are determined within the period of warranty.

The warranty excludes components that are replaced during maintenance such as batteries, etc.

The warranty claim extends to restoring the instrument to readiness for use but not, however, to any further claim for damages. Improper handling or unauthorized opening of the instrument invalidates any warranty claim.

To ascertain the warranty liability, return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

**Copyright**

© Weilheim 2002, WTW GmbH & Co. KG

Reprinting - even as excerpts - is only allowed with the explicit written authorization of WTW GmbH & Co. KG, Weilheim.

Printed in Germany.

<b>1</b>	<b>Overview</b>	<b>5</b>
1.1	Keyboard	6
1.2	Display	7
1.3	Sockets	7
<b>2</b>	<b>Safety</b>	<b>9</b>
2.1	Authorized use	10
2.2	General safety instructions	10
<b>3</b>	<b>Commissioning</b>	<b>13</b>
3.1	Scope of delivery	13
3.2	Initial commissioning	13
<b>4</b>	<b>Operation</b>	<b>15</b>
4.1	Switch on the instrument	15
4.2	Measuring	16
4.2.1	Measuring the pH value	17
4.2.2	Measuring the Redox voltage	18
4.3	Calibrating	19
4.3.1	AutoCal TEC	21
4.3.2	AutoCal DIN	23
4.3.3	ConCal	25
4.4	Reset	28
<b>5</b>	<b>Maintenance, cleaning, disposal</b>	<b>29</b>
5.1	Maintenance	29
5.2	Cleaning	30
5.3	Disposal	30
<b>6</b>	<b>What to do if...</b>	<b>31</b>
<b>7</b>	<b>Technical Data</b>	<b>35</b>
<b>8</b>	<b>Lists</b>	<b>39</b>
<b>9</b>	<b>Appendix</b>	<b>45</b>




# 1 Overview

The compact *inoLab pH Level 1* precision pH meter lets you perform pH measurements rapidly and reliably.

The *inoLab pH Level 1* provides the highest degree of operating comfort, reliability and measuring safety for all applications.

The proven MultiCal<sup>®</sup> calibration procedures and special *AutoRead* function support your work with the pH meter.



- 1 Keypad
- 2 Display
- 3 Sockets




## Note

The measuring instrument can also be delivered as part of a set.


Information on this and other accessories is available in the WTW catalog LABORATORY AND FIELD INSTRUMENTATION or via the Internet.

### 1.1 Keyboard




1	Activate/deactivate AutoRead function
2	Call up calibration procedure
3	Select measuring mode
4	Measuring instrument ON/OFF
5	Reduce values, scroll
6	Increase values, scroll
7	Confirm inputs, start AutoRead

## 1.2 Display



## 1.3 Sockets



Connectors:

1	pH electrode
2	Temperature probe
3	Reference electrode
4	Plug-in power supply (option)



### Caution

Only connect probes to the instrument that cannot feed excessive voltages or currents (> SELV and > circuit with current limiter).

Almost all commercial electrodes - especially WTW electrodes - meet these requirements.





## 2 Safety

This operating manual contains basic instructions that you must follow during the commissioning, operation and maintenance of the pH meter. Consequently, all responsible personnel must read this operating manual before working with the instrument.

The operating manual must always be available within the vicinity of the instrument.

### Target group

This measuring instrument was developed for use in the laboratory.

Thus, we assume that, as a result of their professional training and experience, the operators will know the necessary safety precautions to take when handling chemicals.

### Symbols used



#### Caution

indicates instructions that have to be followed to prevent damage to your instrument.



#### Warning

indicates instructions that have to be followed to protect yourself and the instrument from dangerous electrical voltage.



#### Note

Indicates notes that draw your attention to special features.



#### Note

Indicates cross-references to other documents, e.g. application reports, operating manuals of combination electrodes, etc.

### 2.1 Authorized use

This instrument is authorized exclusively for pH and Redox measurements in the laboratory.

The technical specifications as given in chapter 7 TECHNICAL DATA must be observed. Only the operation and running of the measuring instrument according to the instructions given in this operating manual is authorized.

Any other use is considered **unauthorized**.

### 2.2 General safety instructions

This instrument is constructed and tested in compliance with the EN 61010-1 safety regulations for electronic measuring instruments.

It left the factory in a safe and secure technical condition.

#### Function and operational safety

The smooth functioning and operational safety of the instrument can only be guaranteed if the generally applicable safety measures and the specific safety instructions in this operating manual are followed.

The smooth functioning and operational safety of the instrument can only be guaranteed under the climatic conditions specified in chapter 7 TECHNICAL DATA.

If the instrument was transported from a cold environment to a warm environment, the formation of condensate can lead to the faulty functioning of the instrument. In this event, wait until the temperature of the instrument reaches room temperature before putting the instrument back into operation.



#### Caution

The instrument is only allowed to be opened by personnel authorized by WTW.

**Safe operation**

If safe operation is no longer possible, the instrument must be taken out of service and secured against inadvertent operation.

Safe operation is no longer possible if:

- the instrument has been damaged in transport
- the instrument has been stored under adverse conditions for a lengthy period of time
- the instrument is visibly damaged
- the instrument no longer operates as described in this manual.

If you are in doubt contact the supplier of the instrument.

**Obligations of the operator**

The operator of this measuring instrument must ensure that the following laws and guidelines are observed when using dangerous substances:

- EEC directives for protective labor legislation
- National protective labor legislation
- Safety regulations
- Safety datasheets of the chemical manufacturer.



### 3 Commissioning

#### 3.1 Scope of delivery

- Laboratory measuring instrument, inoLab pH Level 1
- Operating manual and short manual
- 4 x type AA Mignon 1.5 V batteries



#### 3.2 Initial commissioning

Perform the following activities:




- Set the °C or °F
- Set the resolution
- Connect the plug-in power supply (optional).

#### Setting °C or °F





The temperature can be displayed in °C or in °F . In the delivery condition, the measuring instrument is preset to °C. To change the unit proceed as follows:

1	Switch the measuring instrument off.
2	Press and hold down the  key.
3	Press the  key.



4	Toggle between °C and °F by pressing   .
5	Confirm with  . The measuring instrument switches to the measuring mode.

### Setting the resolution

1	Press and hold down the  key.
2	Press the  key. The measured values are displayed with a high resolution, e. g. pH = 4.012.
3	Press the  and keys  again. The measured values are displayed with a low resolution, e. g. pH = 4.01.

### Connecting the plug-in power supply (optional)




#### Warning

The line voltage on site must lie within the input voltage range of the original plug-in power supply unit (see chapter 7 TECHNICAL DATA).



#### Caution


Use original plug-in power supplies only (see chapter 7 TECHNICAL DATA).



1	Insert the plug (1) into the socket (2) of the pH meter.
2	Connect the original WTW plug-in power supply (3) to an easily accessible mains socket.

## 4 Operation

### 4.1 Switch on the instrument

1	Place the instrument on a flat surface and protect it against intense light and heat.
2	Press the  key. The <i>Display test</i> appears briefly on the display. The instrument then switches automatically to the previously selected measuring mode.
3	Connect the pH electrode to the instrument. The measuring instrument is ready for operation.



#### Note

The instrument has an energy saving feature to avoid unnecessary battery depletion.


The energy saving feature switches the instrument off if no key has been pressed for an hour.

The energy saving feature is not active if the instrument is supplied by the plug-in power supply.

## Preparatory activities

### 4.2 Measuring

Perform the following activities when you want to measure:

1	Connect the pH electrode to the instrument.
2	Adjust the temperature of the buffer or test solutions or measure the current temperature if the measurement is made without a temperature probe.
3	Calibrate or check the instrument with the pH electrode.
4	Select the measuring mode by pressing  .



#### Note

Incorrect calibration of the pH electrode will result in incorrect measured values. Therefore, regularly perform calibration before measuring.

## Temperature probe



Measurements can be performed with and without a temperature probe. A connected temperature probe is indicated by TP on the display.



#### Note



The pH meter automatically recognizes the type of the temperature probe used. As a result, you can connect electrodes with the NTC30 or Pt1000.

The temperature measurement is absolutely essential for a reproducible pH measurement. If the measurement is made without a temperature probe, proceed as follows:

1	Determine the current temperature using a thermometer.
2	Set up the temperature by pressing   .





#### Note

When calibrating without a temperature probe, set up the current temperature of the respective buffer solution manually by pressing the   keys.



### 4.2.1 Measuring the pH value

1	Perform the preparatory activities according to section 4.2.
2	Immerse the pH electrode into the test sample.
3	Press the  key until <i>pH</i> appears in the status display. The pH value appears on the display.








#### AutoRead AR (Drift control)

The *AutoRead* function (drift control) checks the stability of the measurement signal. The stability has a considerable effect on the reproducibility of the measured values.

For identical measurement conditions, the following criteria apply:

- pH value: better than 0.02 (setting time: > 30 s)

1	Call up the pH measuring mode by pressing  .
2	Activate the AutoRead function by pressing  . The current measured value is frozen (Hold function).
3	Start the AutoRead function by pressing  . AR flashes on the display until a stable measured value is reached.
4	If necessary, start the next AutoRead measurement by pressing  .
5	To terminate the AutoRead function: Press the  key.




### Note


The current AutoRead measurement (with acceptance of the current value) can be terminated at any time by pressing



### 4.2.2 Measuring the Redox voltage

The instrument can measure the Redox voltage (mV) of a solution when connected with a Redox electrode, e.g. SenTix ORP.

1	Perform the preparatory activities according to section 4.2.
2	Immerse the Redox electrode into the test sample.
3	Press the  key repeatedly until U appears in the status line. The Redox voltage (mV) of the sample appears on the display.
4	Wait for a stable measured value.



### Note

Redox electrodes are not calibrated. However, you can check Redox electrodes using a test solution.

### 4.3 Calibrating

#### Why calibrate?

pH electrodes age. This changes the asymmetry (zero point) and slope of the pH electrode. As a result, an inexact measured value is displayed. Calibration determines the current values of the asymmetry and slope of the electrode and they are stored in the instrument. Thus, you should calibrate at regular intervals.

#### When to calibrate?

- After connecting another electrode
- When the sensor symbol flashes:
  - after expiry of the calibration interval
  - after a voltage interruption, e.g. battery change

You can choose between 3 calibration procedures:

#### AutoCal TEC

is specially adapted to the WTW technical buffer solutions as a fully automatic two-point calibration. The buffer solutions are automatically recognized by the instrument.

#### AutoCal DIN

is specially adapted to permanently programmed buffer solutions according to DIN 19266 as a fully automatic two-point calibration. The buffer solutions are automatically recognized by the instrument.

#### ConCal

is the conventional two-point calibration with 2 freely selectable buffer solutions or single-point calibration as the rapid method.






#### AutoRead

In calibration using AutoCal TEC and AutoCal DIN, the *AutoRead* function is automatically activated. The current AutoRead measurement (with acceptance of the current value) can be terminated at any time by pressing




### Calibration evaluation

After the calibration, the instrument automatically evaluates the current status. The asymmetry and slope are separately evaluated. The worst evaluation appears on the display.

Display	Asymmetry [mV]	Slope [mV/pH]
	-15 ... +15	-60.5 ... -58
	-20 ... +20	-58 ... -57
	-25 ... +25	-61 ... -60.5 or -57 ... -56
  Clean the electrode according to the sensor operating manual	-30 ... +30	-62 ... -61 or -56 ... -50
  Clear the fault according to chapter 6 WHAT TO DO IF...	< -30 or > 30	< -62 or > -50

### Preparatory activities

1	Switch on the instrument by pressing  .
2	Connect the pH electrode to the instrument.
3	Keep the buffer solutions ready.
4	Adjust the temperature of the solutions and measure the current temperature if the measurement is performed without a temperature probe.


### 4.3.1 AutoCal TEC


Use any two of the WTW technical buffer solutions for this procedure.







#### Note

Steps 2 and 6 are not required if you use a temperature probe.






- 1 Press the  key repeatedly until the AutoCal TEC function display appears.



- 2 If necessary, set the temperature of the buffer solution by pressing  .
- 3 Submerge the pH electrode in the first buffer solution.
- 4 Press the  key.  
AR flashes on the display.  
The electrode voltage (mV) appears on the display.  
As soon as a stable value is recognized, Ct2 appears.



- 5 Thoroughly rinse the electrode with distilled water.

6	If necessary, set the temperature of the second buffer solution by pressing   .
7	Submerge the electrode in the second buffer solution.
8	Press the  key. AR flashes on the display. The electrode voltage (mV) appears on the display. As soon as a stable value is recognized, AR disappears. The sensor symbol shows the electrode evaluation after the two-point calibration. The value of the slope (mV/pH) appears on the display.
9	Press the  key. The value of the asymmetry (mV) appears on the display.
10	To return to the measuring mode: Press the  key.


### 4.3.2 AutoCal DIN


Use two different DIN buffer solutions (type A, C, D or F with the pH values 1.679, 4.006, 6.865, 9.180) for this procedure.







#### Note

Steps 2 and 6 are not required if you use a temperature probe.






- 1 Press the  key repeatedly until the AutoCal DIN function display appears.



- 2 If necessary, set the temperature of the buffer solution by pressing  .
- 3 Submerge the pH electrode in the first buffer solution.
- 4 Press the  key.  
AR flashes on the display.  
The electrode voltage (mV) appears on the display.  
As soon as a stable value is recognized, Cd2 appears.



- 5 Thoroughly rinse the electrode with distilled water.

6	If necessary, set the temperature of the second buffer solution by pressing   .
7	Submerge the electrode in the second buffer solution.
8	Press the  key. AR flashes on the display. The electrode voltage (mV) appears on the display. As soon as a stable value is recognized, AR disappears. The sensor symbol shows the electrode evaluation after the two-point calibration. The value of the slope (mV/pH) appears on the display.
9	Press the  key. The value of the asymmetry (mV) appears on the display.
10	To return to the measuring mode: Press the  key.



### Two-point calibration

#### 4.3.3 ConCal


Use two buffer solutions for this procedure:


- pH  $7.0 \pm 0.5$
- any other buffer solution











#### Note







Steps 2 and 9 are not required if you use a temperature probe.

- 1 Press the  key repeatedly until the ConCal function display appears.



- 2 If necessary, set the temperature of the buffer solution by pressing  .
- 3 Submerge the pH electrode in the pH  $7.0 \pm 0.5$  buffer solution.
- 4 Press the  key.  
The measured pH value appears on the display.
- 5 Set the nominal pH value of the buffer solution (at the current temperature) by pressing the   keys.
- 6 Press the  key.  
The value of the asymmetry (mV) and the sensor symbol appear on the display.
- 7 Press the  key.  
SLO(pe) appears on the display.



8	Thoroughly rinse the electrode with distilled water.
9	If necessary, set the temperature of the second buffer solution by pressing   .
10	Submerge the electrode in the second buffer solution.
11	Press the  key. The second measured pH value appears on the display.
12	Set the nominal pH value of the second buffer solution (at the current temperature).
13	Press the  key. The value of the slope (mV/pH) appears on the display. The sensor symbol shows the evaluation of the electrode after the two-point calibration.
14	Press the  key. The value of the asymmetry (mV) appears on the display again.
15	To return to the measuring mode: Press the  key.

## Single-point calibration



Use a buffer solution in the range  $\text{pH} = 7.0 \pm 0.5$  for this procedure.


### Note


Only the electrode asymmetry is determined in single-point calibration. The slope of the last two-point calibration is retained.










### Note

Step 2 is not required if you use a temperature probe. The TP display indicates an active temperature measurement.

- 1 Press the  key repeatedly until the ConCal function display appears.



- 2 Set the temperature of the buffer solution by pressing  .
- 3 Submerge the pH electrode in the buffer solution.
- 4 Press the  key.  
The measured pH value appears on the display.
- 5 Set the nominal pH value of the buffer solution (at the current temperature) by pressing the   keys.
- 6 Press the  key.  
The value of the asymmetry (mV) and the sensor symbol for the evaluation of the electrode appears on the display.
- 7 To return to the measuring mode: Press the  key.



### Basic settings


#### 4.4 Reset




The following functions are reset (initialized) to the values they had on delivery:

Measuring mode	pH
Asymmetry	0 mV
Slope	-59.16 mV/pH
Calibration procedure	AutoCal TEC
Temperature, manual	25 °C
Resolution of pH display	0.01

Proceed as follows:

- 1 Press and hold down the  key.
- 2 Press the  key.




- 3 Toggle between no and yes by pressing  .  
yes: reset parameters.  
no: retain settings.
- 4 Confirm with .  
The instrument changes automatically to the pH measuring mode.

## 5 Maintenance, cleaning, disposal

### 5.1 Maintenance

The measuring instrument is almost maintenance-free. The only maintenance task is replacing the batteries:

1	Open the battery compartment (1) on the underside of the instrument.
2	Remove the four batteries from the battery compartment.
3	Insert four new batteries (Type Mignon AA) into the battery compartment.
4	Close the battery compartment (1).



#### Caution

Make sure that the poles of the batteries are the right way round.

The  $\pm$  signs in the battery compartment must correspond to the  $\pm$  signs on the batteries.

Only use leakproof alkaline manganese batteries.



#### Note

See the relevant operating manual of the electrode for instructions on maintenance.

### 5.2 Cleaning

Occasionally wipe the outside of the measuring instrument with a damp, lint-free cloth. Disinfect the housing with isopropanol as required.



#### Caution

The housing is made of synthetic material (ABS). Thus, avoid contact with acetone or similar detergents that contain solvents. Remove any splashes immediately.

### 5.3 Disposal

#### Packing

The measuring instrument is sent out in a protective transport packing.

We recommend: Keep the packing material. The original packing protects the instrument against damage during transport.

#### Batteries

This note refers to the battery regulation that applies in the Federal Republic of Germany. We would ask end-consumers in other countries to follow their local statutory provisions.



#### Note

In compliance with §14 of the BATTERY REGULATION, we would like to point out that this instrument contains batteries. Batteries that have been removed must only be disposed of at the recycling facility set up for this purpose or via the retail outlet.

It is illegal to dispose of them in household refuse.

#### Measuring instrument

Dispose of the measuring instrument as electronic waste at an appropriate collection point. It is illegal to dispose of them in household refuse.

## 6 What to do if...

Error message,  
OFL

Cause	Remedy
pH electrode:	
– Not connected	– Connect electrode
– Air bubbles in front of the diaphragm	– Remove air bubbles
– Air in the diaphragm	– Extract air or moisten diaphragm
– Cable broken	– Replace electrode
– Gel electrolyte dried out	– Replace electrode

Error message,  
E3

Cause	Remedy
pH electrode:	
– Diaphragm contaminated	– Clean diaphragm
– Membrane contaminated	– Clean membrane
– Moisture in the plug	– Dry plug
– Electrolyte obsolete	– Replenish electrolyte or replace electrode
– Electrode obsolete	– Replace electrode
– Electrode broken	– Replace electrode
Measuring instrument:	
– Incorrect calibration procedure	– Select correct procedure
– Incorrect solution temperature (without temperature probe)	– Set up correct temperature
– Socket damp	– Dry socket

## What to do if...


<b>No stable measured value</b>	Buffer solutions:	
	– Incorrect buffer solutions	– Change calibration procedure
	– Buffer solutions too old	– Only use once. Note the shelf life
	– Buffer solutions depleted	– Change solutions
	<b>Cause</b>	<b>Remedy</b>
	pH electrode:	
	– Diaphragm contaminated	– Clean diaphragm
	– Membrane contaminated	– Clean membrane
	Sample:	
	– pH value not stable	– Measure with air excluded if necessary
	– Temperature not stable	– Adjust temperature if necessary
<b>LoBat</b>	Electrode + sample:	
	– Conductivity too low	– Use suitable electrode
	– Temperature too high	– Use suitable electrode
	– Organic liquids	– Use suitable electrode
	<b>Cause</b>	<b>Remedy</b>
	– Batteries almost depleted	– Replace batteries (see section 5.1 MAINTENANCE)




**Obviously incorrect  
measured values**

Cause	Remedy
pH electrode:	
– pH electrode unsuitable	– Use suitable electrode
– Temperature difference between buffer and sample too large	– Adjust temperature of buffers or samples
– Measuring procedure not suitable	– Follow special procedure

**Instrument does not  
react to keystroke**

Cause	Remedy
– Operating state undefined or EMC electric stress unallowed	– Processor reset: Press the  key and switch on instrument

**You would like to  
know which soft-  
ware version is in  
the instrument**

Cause	Remedy
– e.g. question of the WTW service department	– Press the  key and switch on instrument. The software version is displayed.



## 7 Technical Data

<b>Ambient temperature</b>	Storage temperature	- 25 °C ... + 65 °C
	Operating temperature	0 °C ... + 55 °C
	Allowable relative humidity	Annual mean: < 75 % 30 days/year: 95 % Other days: 85 %
<b>Measuring ranges and resolution</b>	pH	- 2.000 ... + 19.999 - 2.00 ... + 19.99
	U [mV]	- 999.9 ... + 999.9 - 1999 ... + 1999
	T [°C]	- 5.0 ... + 105.0
	T [°F]	+ 23.0 ... + 221.0
<b>Accuracy (± 1 digit)</b>	pH (in the measuring range of 2 pH units around the calibration point)	± 0.005 (at operating temperature + 15 °C ... + 35 °C ) ± 0.01
	U [mV]	± 0.3 (at + 15 °C ... + 35 °C) ± 1
	T [°C]	NTC 30: ± 0.1 PT 1000: ± 0.5 at 0 °C ... 15 °C ± 0.1 at 15 °C ... 35 °C ± 1 at 35 °C ... 55 °C
	T [°F]	NTC 30: ± 0.2 PT 1000: ± 0.9 at 32 °F ... 59 °F ± 0.2 at 59 °F ... 95 °F ± 1.8 at 95 °F ... 131 °F

## Technical Data

---

### Dimensions and weight

Length [mm]	230
Width [mm]	210
Height [mm]	70
Weight [kg]	Approx. 0.850

### Energy supply

Batteries	4 x 1.5 V AA type alkaline manganese batteries
Runtime	Approx. 3000 operating hours
Mains power supply (option)	<p>Connection max. overvoltage category II (valid for all plug-in power supplies):</p> <p>Plug-in power supply (Euro plug): FRIWO FW1199, 11.7864 Friwo Part. No. 1762613 Input: 230 V ~ / 50 Hz / 5.6 VA Output: 12 V = / 130 mA / 1.56 VA</p> <p>Plug-in power supply (US plug): FRIWO FW1199, 11.7880 Friwo Part. No. 1794043 Input: 120 V ~ / 60 Hz / 6 VA Output: 12 V = / 150 mA</p> <p>Plug-in power supply (UK plug): FRIWO FW3288, 11.8453 Friwo Part No. 1816491 Input: 230V ~ / 50 Hz / 23 VA Output: 12 V = / 130 mA / 1,56 VA</p>

**Guidelines  
and norms used**

EMC

E.C. guideline 89/336/EEC  
EN 61326-1:1997  
EN 61000-3-2 A14:2000  
EN 61000-3-3:1995  
FCC Class A

Instrument safety

E.C. guideline 73/23/EEC

Protective class

3, EN 61010-1 A2:1995

Climatic class

2, VDI/VDE 3540

**Test marks**

TÜV GS, UL/CUL, CE



## 8 Lists

This chapter provides additional information and orientation aids.

### **Abbreviations**

The list of abbreviations explains abbreviations that appear on the display or when dealing with the instrument.

### **Specialist terms**

The glossary briefly explains the meaning of the specialist terms. However, terms that should already be familiar to the target group are not described here.

### **Index**

The index helps you find the topics that you are looking for.

**Abbreviations**

AR	AutoRead (drift control)
ARng	Automatic range switching Measuring instrument measures with highest resolution
ASY	Asymmetry
AutoCal DIN	Automatic calibration with DIN buffer solutions
AutoCal TEC	Automatic calibration with WTW technical buffer solutions
Cal	Calibration
Cd...	Calibration with DIN buffer solutions (acc. to DIN 19 266)
ConCal	Conventional one/two point calibration
Ct...	Calibration with WTW technical buffer solutions
E3	Error message (see WHAT TO DO IF ...)
InI	Initialization Resets individual basic functions to the status they had on delivery
LoBat	Low Battery Batteries are almost empty
mV	Voltage unit
mV/pH	Unit of the electrode slope
OFL	Overflow Display range exceeded
pH	pH value
S	Slope
SELV	Safety Extra Low Voltage



SLO	Slope Slope setting on calibration
TP	Temperature probe Temperature measurement active
U <sub>ASY</sub>	Asymmetry potential
°C	Temperature unit, °Celsius
°F	Temperature unit, Fahrenheit

## Glossary

<b>Asymmetry</b>	Zero point of a pH electrode.
<b>AutoRead</b>	Monitors the electrode drift and releases the measured value only after the stability criterion has been reached. In this way, this procedure ensures the highest degree of precision and reproducibility.
<b>Buffer solution</b>	Stable solution with a precisely known pH value.
<b>Diaphragm</b>	Contact point between the reference electrolytic solution and the sample.
<b>Drift control</b>	See AUTOREAD.
<b>MultiCal<sup>®</sup></b>	Group term for the various WTW calibration procedures used for automatic calibration in buffer solutions.
<b>Redox voltage</b>	Potentiometric quantity.
<b>Resolution</b>	Number of decimal places that appear for a measured value.
<b>Slope</b>	Specifies the voltage change per pH unit.
<b>Test solution</b>	Stable solution with a precisely known Redox voltage.

**Index**
**A**

asymmetry 19  
 authorized use 10  
 AutoCal DIN 19, 23  
 AutoCal TEC 19, 21  
 AutoRead 17

**B**

basic settings 28  
 battery compartment 29

**C**

calibrating 19  
 calibration evaluation 20  
 calibration procedures 19  
 ConCal 19, 25

**D**

delivery condition 28  
 display 7  
 drift control 17

**E**

electrode symbol 7  
 energy saving feature 15  
 error messages 31

**I**

initial commissioning 13  
 initialize 28

**K**

keys 6

**L**

LoBat 32

**O**

operational safety 10

**P**

place of the instrument 15  
 plug-in power supply 14

**R**

Redox electrode 18  
 Redox voltage 18  
 replacing the batteries 29  
 reset 28

**S**



safety 9  
 safety precautions 9  
 scope of delivery 13  
 setting the resolution 14  
 setting the temperature unit 13  
 single-point calibration 19  
     ConCal 27  
 slope 19  
 sockets 7

**T**

temperature probe 16  
 two-point calibration 19  
     AutoCal DIN 24  
     AutoCal TEC 22  
     ConCal 25



## 9 Appendix

<b>SENTON</b>	
<b>CERTIFICATE</b>	
<b>to the Electromagnetic Compatibility</b>	
<b>to Test Report No. 52501-90109/-1</b>	
<b>EUT:</b>	InoLab pH Level 2 with printer
<b>Uniform EMC-design:</b>	InoLab pH Level 1, InoLab pH Level 2, InoLab pH Level 2 with printer, InoLab pH/ION Level 2, InoLab pH/ION Level 2 with printer, CG842, CG843, CG843P
<b>Applicant:</b>	WTW GmbH
<b>Regulations:</b>	EN 50081-1:1992 EN 50082-1:1992 FCC part 15 subpart B
<b>Test result:</b>	
The tested sample is in compliance with the RFI requirements and the immunity requirements according to above referenced regulations.	
The following severity levels have been achieved:	
RFI-Emissions:	Requirements according to EN 50081-1:1992 FCC part 15 subpart B limit class A
Immunity:	Electrostatic Discharge according to IEC 801-2:1984 Air discharge: 8 kV  Radiated Electromagnetic Fields according to IEC 801-3:1984 27 - 500 MHz: 3 V/m  Electrical Fast Transients (Burst) according to IEC 801-4:1988 Power line: 1 kV; data line: 0.5 kV
 Roidt Johann Senton GmbH	March 10, 1999 Date
	DAR-Registration No. TTI-P-G 062/94-30 DAR-Registration No. TTI-P-G 109/95-20
SENTON GmbH - EMV-Prüfzentrum - Äußere Frühlingsstraße 45 - D-94315 Straubing - Tel. 09421/55 22-0	



### **Note**

If you need further information or application notes you can ask for:

- Application reports
- Primers
- Safety data sheets.

You can find information on available literature in the WTW catalog, LABORATORY AND FIELD INSTRUMENTATION or via the Internet.